
# Simbiosi -Pre-Milled®Prodotti per protesi integrale Products for integral prosthesis

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Simbiosi Pre-Milled® can be used for the construction of custom abutments, suitable for each patient, while treating final dental restorations or temporary teeth.

In the pre-milling technique of Simbiosi Pre-Milled®, the user builds the morphology of the abutment with any CAD system designer, neglecting the realization of each individual implant's engagement. The system involves the construction of the abutment's morphology milling a Pre-Milled® that is already provided with a pre-processed engagement, to always guarantee a perfect fit.

The product Simbiosi Pre-Milled® allows you to always have custom trans-mucosal routes without any effort, regardless of the implant position. Therefore it allows time savings in producing abutments and lets you orient every insertion axis of the secondary prosthesis. To realize your abutment with Simbiosi Pre-Milled® becomes simple, upgraded, and cheap.

Simbiosi Pre-Milled® is available in different sizes, diameters and shapes, respectively compatible with the implants' manifacturers, with special diameter and particular implant system.

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

Simbiosi produces three types of Scan Body:

Scan Body IPF - Scan Body UNICO - Scan Body TOOTHLESS  $^{\tiny{\textcircled{\tiny 0}}}$ 

With the aid of the scanner (desktop or intraoral), these tools are suitable devices for the detection of the exact position of the endosseus implant that the professional inserted into the patient's mouth. The aim is to physically or virtually reconstruct the exact reproduction of the oral cavity that has to receive the dental prosthesis. The Scan Body is a necessary device and it is solely used for the implant position's detection to build abutment in electronic format using CAD softwares.

Simbiosi Scan Body is available in different sizes, diameters and shapes, respectively compatible with the implants' manifacturers, with special diameter and particular implant system.

## Come e dove usare Simbiosi Pre-Milled How and where to use Simbiosi Pre-Milled

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Pre-milled® system was born from an idea by Simbiosi to satisfy all the needs of the dental prosthesis'construction; the prescription of cemented prosthesis with abutment and subsequent dental element, is commonly used, but hard to find in today's market at competitive prices, or the technician must comply with standard shapes that are difficult to be adapted to the patient. For this reason, Simbiosi has devised Pre-milled® system, directly usable in laboratories with milling machines for metal. By partnering with manufacturers of milling machines, Simbiosi offers the possibility to install your Pre-milled® system in simplicity, certifying the perfect operation. For the correct use of the Pre-milled® you must install the "bracket" support to receive it and then work it with a milling undercut procedure, calculated by the most famous CAM that are currently on the market, combined with milling machines.

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THE INTERNATIONAL CERTIFICATION NETWORK

## **CERTIFICATE**

IQNet and its partner CISQ/CERTIQUALITY S.r.I.

SIMBIOSI S.R.L.

IT - 50053 EMPOLI (FI) - VIA PAOLO VERONESE 2/A
has implemented and maintains a

Quality Management System
withs fulfair to requirements of the following standard
ISO 13485:2003

Code AT 1d
sign, production and sale of medical devices for the dentistry.

in the following operative units
IT - 50035 BMPOLI (Fi) - VNA PLAVE 3/5
IT - CORRIDONIA (MC) - VNA PLAVE 3/5
IT - CORRIDONIA (MC) - VNA PLAVE 3/5
Issued on: 2015-01-08
Certified since: 2015-01-08
Expire on: 2018-01-07
Registration number: IT-995 IT-99561

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